SPECIMEN SUBMISSION GUIDELINES

HIV Diagnostics and Reference Laboratory, Diagnostics and Countermeasures Branch Center for Infectious Disease Research Walter Reed Army Institute of Research 508 Research Dr., Silver Spring, MD

Hepatitis B Surface Antigen Test Request Form

Test Requested	Specimen Requirement		Draw Tube	Shipping Conditions (Check one)		
☐ HBsAg Confirmatory	☐ 3 ml serum (SST Tubes) or ☐ 3 ml plasma (EDTA, Na Heparin, Na Citrate, CPDA and ACD-1 plasma is acceptable.)	to clot for hrs) post-c 10 minutes a swing bu NOTE: To to clot for PPT Tube Stable at r 6 hrs. Cen rotor centra minimur	es – Invert 5X and allow 30 min (no more than 2 collection. Centrifuge s at 1000-1300 RCF in acket centrifuge. ubes MUST be allowed 30 minutes. es – Invert 8-10X. coom temperature up to trifuge in swing-out rifuge at 1100 RCF for m of 10 min. Freeze quot at -20°C.	□ Ambient 15-30°C – SST tube must be received at HDRL within 7 days of collection. □ Refrigerated 2-8°C – SST tube must be shipped in cold box with ice packs and received at HDRL within 7 days of collection. □ Frozen -20°C – Ship frozen aliquoted plasma with dry ice if specimen will be received at HDRL after 7 days of collection.		
Please fill the request form completely to ensure timely specimen processing.						
PATIENT IDENTIFICATION			CONTACT INFORMATION			

PATIENT IDENTIFICATION

Patient identifiers MUST INCLUDE:

Full Name

OR

FMP/SSN

OR

Telephone Number

Fax Number

(Commercial # only; please include area/country code)

Alternate POC Name

Alternate POC Name

PROCESSING LABORATORY (For HDRL use only)					
BARCODE	DATE RECEIVED	QUANTITY & TYPE RECEIVED / INITIALS			

Ship Date:

Alternate POC Phone

Fax/Email a FedEx tracking and/or invoice number to ensure all shipments sent to the HIV Diagnostics and Reference Laboratory are received, IAW CAP GEN.40530
Shipping Address: 9100 Brookville Road, BLDG 508, Silver Spring, MD 20910

Form # TR HBsAg Version: June 2023