SPECIMEN SUBMISSION GUIDELINES

HIV Diagnostics and Reference Laboratory, Diagnostics and Countermeasures Branch Center for Infectious Disease Research Walter Reed Army Institute of Research 508 Research Dr., Silver Spring, MD

Molecular Aptima HIV-1 and HCV Qualitative/Quantitative Dx Test Request Form

Test Requested	Specimen Requirement	Draw Tube	Shipping Conditions (Check one)	
Qual / Quant HIV-1 (Viral Load) Qual / Quant HCV (Viral Load)	☐ 3 ml plasma (PPT or EDTA)	PPT tubes – Invert 8-10X. Spin tubes within 2 hrs of collection. Centrifuge ≤1300 RCF for 10 minutes. Freeze plasma aliquot at less than -20°C. EDTA tubes – Store blood at 25°C until centrifuged. Centrifuge ≤1300 RCF for 10 minutes within 4 hours of blood collection. Freeze plasma aliquot at less than -20°C.	☐ Frozen -20°C — Ship frozen aliquoted specimen with dry ice if specimen will be received at HDRL after 24 hours of collection. RECOMMENDED: Ship on 6 lbs. dry ice in case of shipment delay.	
Please fill the request form completely to ensure timely specimen processing.				

rease in the request form completely to ensure timely specimen processing.				
PATIENT IDENTIFICATION	CONTACT INFORMATION			
Patient identifiers <u>MUST INCLUDE</u> :	POC			
Full Name	Physician Name			
DoDIDOR	Center Address			
FMP/SSN				
Specimen Draw Date / Time: Ship Date:	Fax Number (Commercial # only; please include area/country code)			
PROCESSING LABORATORY (For HDRL use only)				
BARCODE	DATE RECEIVED QUANTITY & TYPE RECEIVED / INITIALS			

Fax/Email a FedEx tracking and/or invoice number to ensure all shipments sent to the HIV Diagnostics and Reference Laboratory are received, IAW CAP GEN.40530
Shipping Address: 9100 Brookville Road, BLDG 508, Silver Spring, MD 20910

Form # TR Aptima Version: June 2023